## Temp S

Handmixable Powder/Liquid System for Temporary Creating Crowns and Bridges

**Temp S** is a special handmixable two component powder/liquid system for creating temporary crowns or bridges. As a temporary crown or bridge it protects the prepared teeth against external influences and preserves the occlusion. It shows a good chewing stability and fits perfect.

### Application

#### Impression taking

**Before preparing of stumps or before a planned extraction of a tooth** a situation impression is made by addition curing silicones. Interdental areas are carved out, in molar areas with teeth absent it is may be necessary to cut a groove in the impression between the abutments to create a bridge-like connection between the tooth units.

#### Mixing

Poor into the mixing pot a quantity of liquid, which is adequate to the volume of the temporary, which has to be prepared. Add the powder just up a middle flowing consistency results (mixing ratio powder/liquid = 2:1).

#### Application

Apply **Temp S** into the situation impression.

#### Forming of the temporary crowns or bridges by indirect technique

The situation impression - used as a tray and loaded with **Temp S** like in a putty/wash technique - is seated onto the prepared areas of the model. After about **6-7 minutes** at room temperature the material shows a hardened but still elastic condition and can be removed from the model.

#### Post curing and working out

After removal of the temporary crown or bridge from the situation impression and tooth model excess material and proximal undercuts are removed. Preferably, the temporary is then post cured in hot water (e.g. in a hot cure polymerisation device). If this is not possible or not desired, after **13 minutes**, the temporary acrylic can be worked out with rotary instruments and can be polished to high gloss like a microfill or a small particle hybrid composite material.

#### Note:

The oxygen inhibited smear layer on the surface, caused by the acrylic system, can not be fully avoided and should be removed before working out. It can easily be removed by soft paper tips soaked with alcohol or other suitable solvents.

#### Cementing of the temporary

**Temp S** temporary crowns or bridges preferably should be cemented with an eugenol-free temporary cement. In case of using eugenol-containing cements it should be considered that later used acrylate based composite cements could be hindered in hardening (inhibition of curing reaction by traces of remained eugenol).

#### Repairs

**Temp S** temporaries show high mechanical strength. However, if a temporary breaks the following procedures are recommended:

##### a) Fracture shortly after production:

Both ends of the fracture are refixed with freshly mixed **Temp S**.

##### b) Fracture of longer existing temporary:

The areas of fracture are cleaned and roughened and provided with some mechanical retentions. Thus prepared, the fractured areas can be joint with freshly mixed **Temp S**. For enhancing the curing, place for some minutes in warm water (50 °C).

#### Special hints

Non hardened **Temp S** can easily removed with alcohol, acetone or other suitable solvents.

### Precautions

In some patients **Temp S** may cause sensibilization. In these cases **Temp S** should not be used any more.

#### Warnings

The **Temp S** contains methacrylates. **Temp S** Liquid is irritating. Contact of the material with skin, mucous membrane or eyes should be avoided. After contact with skin wash off with plenty of water. After eye contact rinse out with plenty of water and see a physician.

#### Storage

Do not store above 25° C (77 °F)! **Temp S** shall not be used after expiry date!

#### Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico's sole obligation shall be repair or replacement of the Bisico product.

#### Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

#### Time Table of Application and Curing of Temp S

START ↓ 1 : 30 min ↓ 6:00 - 7:00 min ↓ 13:00 min.	Start of mixing: mixing ratio: powder/liquid 2:1  Working time (from begin of mixing)  Setting time (23°C) on model to elastic state  Earliest start for trimming and polishing (full hardened state)
---	---

#### Keep away from children! For dental use only!

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