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Repo-Core

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Instructions for Use

Repo-Core DC

Dual Curing Flowable Microhybrid Composite for Core Build-ups and Cementation of Posts

Product description

Repo-Core DC is a dual curing radiopaque flowable microhybrid composite based material for core-build-up and cementation of posts. Repo-Core DC is also suitable for cementing of crowns and bridges, inlays and onlays. For all indications the use of a dual cure bonding agent (eg. E-Bond DC) is required before application of the composite.

Repo-Core DC is based on poly- and difunctional methacrylates and inorganic filler particles of 0.05-2.5 µm. The total filler content is 66 % by weight and 46 % by volume. Delivered in auto-mixing 1:1 syringes it can be easily dispensed and applied directly. Repo-Core DC exhibits a short setting time without high heat generation. The dual cure properties enables the dentist, also to carry out cementations and core build ups in cases where a light cure cannot be guaranteed to be sufficient. Repo-Core DC meets the requirements of DIN EN ISO 4049, type 2, class 3.

Indications/Intended use

- Core build-ups
- Cementing of posts
- Cementing of crowns and bridges, inlays and onlays

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

The placement of Repo-Core DC is contraindicated if a dry working area or the recommended application technique are not possible.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparing the Syringe

Remove the cap of the syringe and throw it away (do not use it again!). It is replaced by a supplied 1:1 mixing cannula. Turn the cannula 90° until it locks in position. The material is now ready for application.

Note:

Store used syringe with fixed used mixing cannula. Discard the first 2-3 mm of the extruded material. This has to be done for each new mix.

Store used syringe with fixed used mixing cannula in the dark.

The working time (23°C (74°F)) of Repo-Core in the self cure mode is 1:30 minutes from start of mixing.

Application

1. Post Cementation

1.1. Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

1.2. Root Canal Preparation

Refer to directions of the selected post manufacturer.

Before starting the preparation clean the tooth from residues. Prepare and clean the root canal with e.g. sodium hypochlorite solution, rinse and remove excess solution from the canal with a soft paper tip.

Recommended bonding agents for conditioning are:

- a suitable dual cure bonding system for the Total-Etch-Technique in combination with Bisico Best Etch
- E-Bond DC, a self-etching dual cure adhesive

Apply the bonding agent according to the corresponding instructions.

1.3. Post Cementation

Prepare the selected post according to manufacturer directions.

Repo-Core DC is applied into the prepared root canal and onto the post. Seat the post careful into the canal and maintain firm pressure until the post is seated.

Repo-Core DC self-cures within 3:30 minutes. For post stabilization light cure the coronal part of the cemented post for 20 seconds with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm2.

As soon as the Repo-Core DC has set proceed with the core-build-up procedure.

2. Core-Build-Up

2.1 Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

2.2. Cavity Preparation

Remove all existing old restorations and decay from the tooth. If necessary place any pins or posts. Refer to directions of the selected post manufacturer.

2.3. Pulp Protection

For pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material (eg. Bisico Bond).

2.4. Application of a Bonding Agent

Recommended bonding agents for conditioning are:

- a suitable dual cure bonding system for the Total-Etch-Technique in combination with Bisico Best Etch
- E-Bond DC, a self-etching dual cure adhesive

Apply the bonding agent according to the corresponding instructions.

Note:

It is essential that the primed dentine and enamel surfaces are dry and contaminant free for the application of Repo-Core DC.

2.5. Application of Repo-Core DC

Place the mixing cannula directly into the preparation and press out the paste. Repo-Core DC is automatically mixed when dispensed with slight and even pressure. Filling should occur from bottom upwards to prevent air voids. To facilitate placement of Repo-Core DC place a matrix band around the prepared tooth.

Repo-Core DC may be contoured by using a suitable modeling instrument.

Place Repo-Core DC directly into the preparation and allow the system to self cure for 3:30 minutes. After that the material should be light-cured (40 seconds) with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm2. With this technique an optimum of physical properties will be obtained.

An explorer can be used to test that Repo-Core DC has completely set. Remove the matrix not earlier than the material has set.

Final core preparation on Repo-Core DC can be carried out by using crown preparation bur.

Additional Notes/Warnings

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of the composite.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Color stability meets the requirements of DIN EN ISO 4049.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

Composition

Dimethacrylates, dental glass, silicone dioxide, photo initiators, catalysts

Storage

Do not store above 20 °C (68 °F). Protect from direct sunlight.Store unopened material in the refrigerator. Opened cartridges have to be used up within 3 months.

Do not use after expiry date.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico GmbH's sole obligation shall be repair or replacement of the Bisico GmbH product.

Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.