Micro Esthetic

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Micro Esthetic

Light curing Microhybrid composite

Micro Esthetic is a visible light cure Microhybrid composite for anterior and posterior restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0.05-1.5 µm. The total filler load is 81 %, the total filler volume is 65 %. Radiopacity of Micro Esthetic is 200 % aluminium (radiopacity of enamel:-200 % aluminium, of dentine:- 100 % aluminium).

The composite material meets the requirements of: DIN EN ISO 4049, type 1, class 2, group 1.

Indications

Suitable for all cavities.

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of hard setting calcium hydroxide material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparation of Tooth and Color Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching. Ascertain the tooth shade while teeth are still moist and select the appropriate Micro Esthetic shade.

Micro Esthetic is shaded according to VITA®-Shades.

Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a thin layer of Bisical or a calcium hydroxide material. Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with Best Etch and subsequent application of a light cure bonding agent (e.g. Bond LC) is recommended. Otherwise use a self etching bonding agent (e.g. E-Bond LC) before application of the Micro Esthetic.

Application and Curing of Micro Esthetic

Apply Micro Esthetic in the selected shade. Micro Esthetic is easy to model. Place it by using a non-metallic instrument. Transparent matrix strips may be used.

In case of deep and extended cavities best results are obtained with a light cure flow able composite e.g. Micro Esthetic flow-viscous or Bisico Bulk Fill, followed by the application of Micro Esthetic.

In case of small cavities Micro Esthetic can be applied directly.

For an optimum result apply Micro Esthetic in layers of max. 2 mm.

Light cure each increment separately.

By using a polymerization (wavelength range 400-500 nm) unit with a light intensity of at least 1000mW/cm², cure each increment as follows:

Lighter shades (e.g. A1, A2, B1; B2; clear; bleach; bleach light) 20 sec.
darker shades (e.g. A3.5, A4; B3) 30 sec.

Instructions for use EN

Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

Storage

Do not store above 25°C (78°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light. Capsules are for single use only.

Additional Notes

- The ambient light of the dental lamp may start polymerization of the composite.

- Do not use any resin to adjust viscosity of composite restorative material.

- Avoid contact with skin, mucous membrane and eyes.

- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.

- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, our exclusively remedy and Bisico's sole obligation shall be repair or replacement of the Bisico product.

Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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bisico Bielefelder Dentalsilicone GmbH & Co. KG • Johanneswerkstraße 3 • D-33611 Bielefeld