**Bisico® Bulk Fill****Zweckbestimmung**

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0297

**Bisico<sup>®</sup> Bulk Fill****Intended Purpose**

Dental filling materials on composite basis are used to build up or to preserve tooth structures

**Product description**

**Bisico<sup>®</sup> Bulk Fill** is a light curing, flowable, highly radiopaque (210% Al) composite for posterior restorations, for direct filling and the „bulk fill“ technique.

Bulk Fill Composite is available in syringes and compules. The compules are for single use. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation otherwise cannot be excluded.

**Patient target group**

All patients requiring tooth substance to be replaced or built up in order to restore or maintain dental function.

**Users**

Dental filling materials on composite basis are used in dental practice or in a dental hospital by dental professionals.

**Composition**

Glass powder, silicon dioxide, 1,4-butanediol dimethacrylate

Filler content: 77% by weight (57% by volume) inorganic fillers (0.005 - 40 µm)

**Indications:**

- Fillings with layer thicknesses up to 4 mm in Black's class I
- Base fillings/liners with layer thicknesses up to 4 mm in Black's class II cavities. For the following layers in a class II cavity, only use composites which are suited to filling Black's class II cavities.
- Cavity lining - as a first (bottom) layer in Black's class I and II cavities

**Contraindications**

If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist.

**Warnings**

Contains :di-urethane dimethacrylate, tetramethylene dimethacrylate

Important: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects.

**Safety instructions**

Wear protective gloves. Avoid breathing vapours und spray.

**Interactions with other materials**

Phenolic substances (such as eugenol) inhibit polymerisation. Therefore, do not use any lining materials (such as zinc oxide-eugenol cements) which contain substances of this type.

The dentist should consider known interactions and cross-reactions of the medical product with other materials already in the patient's mouth before using the product.

**1 Application:**

Before beginning the treatment, clean the hard tooth substance with non-fluoride polishing paste.

**1.1. Cavity preparation**

Prepare the cavity with a minimally invasive method as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. Isolate the cavity to keep it dry. It is advisable to place a rubber dam.

**1.2 Pulp protection / liner**

Lining is not necessary if an enamel-dentin adhesive is used. If there are very deep cavities near the pulp, cover corresponding areas with a calcium hydroxide (e.g. Bisical) preparation.

**1.3. Approximal contact areas**

For cavities with approximal sections, place a transparent matrix and fix it in place.

**1.4. Adhesive system**

Etch and bond according to manufacturer's instructions.

**1.5.1. Application of Bisico<sup>®</sup> Bulk Fill (syringes)**

Lock the application tips on the syringe by turning it clockwise. Place layers of max. 4 mm of **Bisico Bulk Fill** directly into the cavity using the curved application tips supplied with the material. For hygiene reasons, these application tips must only be used once and are then discarded. Avoid introducing air bubbles during application with the tip. Ensure the prepared hard tooth tissue is thoroughly wetted.

**1.5.2 Application of Bisico<sup>®</sup> Bulk Fill (compules)**

Insert the compule into the dispenser. Remove the sealing cap. Place the compule in the correct angle towards the cavity. Inject the material into the cavity. Apply slow and steady pressure to the compule. Do not use excessive force! (max. layer thickness 4 mm) To remove the compule from the dispenser after use, retract the plunger. Next, remove the compule.

Please note: For hygiene reasons, compules are intended for single use only.

**1.6. Curing**

Cure each layer after coating for 20 seconds using a commercially available polymerisation unit. The light guide must be held as close as possible to the surface of the filling. A dispersion layer forms on the surface during polymerisation which should not be touched or removed. It bonds to the next composite layer (e.g. a universal or posterior dental composite) or the aesthetic cusp onlay. Alternatively, a thin layer of Bisico Bulk Fill can be applied to exposed dentine as a liner. Polymerisation can be applied incrementally up to 4 mm.

**1.7. Finishing**

**Bisico<sup>®</sup> Bulk Fill** can be finished and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishers, and polishing brushes can be used for finishing. Check the occlusion and articulation and grind so that no early contacts or undesirable paths of articulation remain on the surface of the restoration.

**2. Special information**

- In case of time-consuming restorations, the surgical lamp should be temporarily moved away from the working area to prevent the composite from curing too early or the material should be covered with an opaque foil.
- The curved delivery tips supplied with the kit are intended for single use only since contamination of the materials and spread of germs cannot be ruled out.
- A light polymerisation unit with an emission spectrum in the range of 350 - 500 nm is to be used for polymerisation. The physical properties required are achieved only with properly working lamps. For this reason, it is necessary to regularly check the light intensity according to the manufacturer's instructions.

Light intensity for curing	≥ 1200 mW/cm <sup>2</sup>
Wavelength for curing	350 - 500 nm
Curing time	20 sec.

**Troubleshooting**

Trouble	Cause	Remedy
composite does not cure	luminous intensity of the polymerization unit insufficient	check luminous intensity; replace light source, if necessary
	emitted spectral range of the polymerization unit insufficient	consult manufacturer of polymerization unit; recommended spectral range: 350 - 500 nm
composite seems to be too hard / firm inside the syringe / capsules	material was stored at temperatures below 10 °C for a longer period of time	let composite reach room temperature before use
	syringe was not closed tightly which caused part of the material to cure	close syringe tightly after each use
composite does not cure sufficiently	layer thickness per polymerization cycle too high	keep to max. layer thickness of 4 mm
restoration seems too yellow when compared to color reference	insufficient polymerization of the composite layers	repeat polymerization cycle several times, for a minimum of 20 sec

**Information on storage and handling**

Store at 10 - 25°C (50 - 77°F). Close syringes tightly immediately after use. The material should be at room temperature before use. Retract the syringe spindle slightly after use to prevent the discharge opening from becoming clogged.

**Shelf life**

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

**Side effects**

With proper preparation and use of this medical product, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case

**Disposal**

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

Bisico Bielefelder Dentalsilicone GmbH & Co. KG  
Johanneswerkstraße 3  
D-33611 Bielefeld  
Germany

MANI MEDICAL GERMANY GmbH  
Hertha-Sponer-Straße 2  
61191 Rosbach v.d.Höhe  
Germany

0297