# *bisico* Softbase

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## Softbase

Permanently soft denture relining material based on addition-curing silicone

**Softbase** is a soft liner used to quickly and easily produce soft denture relinings. Softbase helps meet the demanding requirements for enhanced denture wearing comfort. It is a complete system comprising a primer (Softbase Primer) to ensure adhesion to the base of the denture, the actual relining material (Softbase Cartridge) and a 2-component-based gloss seal (Softbase Varnish).

**Softbase** is a permanently soft direct relining material based on addition-curing silicone in cartridges. It is suitable for all PMMA plastics. The relining can be performed intraoral, directly on the patient, as well as extra orally in the laboratory. It meets the requirements of DIN EN ISO 10139-2.

Softbase prevents the formation of pressure points in the area of the upper and lower jaw. It is used to dam the A-line in cases of difficult adhesion and flabby ridges. It forms an elastic margin on upper dentures and promotes positional stability of the dentures through increased adhesion. Sharp-edged alveolar processes are elastically covered.

Softbase can also be used for damming implants and as overdentures during the healing process.

The **Primer** is used for **adhesion** between the silicone and the denture, and the **Softbase Varnish** is used for **smoothing roughened ground surfaces** of the soft relining material.

### Contraindications

In the case of sensitivity to one or more components of this product, further use of the product should be discontinued.

### Directions for preparing Softbase

#### 1. Preparing the denture

First create the necessary free space for the denture relining. To do this, reduce the base surface by 1-2 mm by grinding or milling and remove any old relining that may be present. In addition, at the vestibular, lingual and dorsal edges, incorporate a step with a depth of at least 1 mm at a height of approx. 2-3 mm. The edges of the denture must be rounded and the surfaces to be relined **must be roughened**.

Following pre-treatment, thoroughly clean the denture and remove traces of grease using pure surgical alcohol.

#### 2. Applying Softbase Primer for adhesion

Coat all areas of the cleaned, ground surfaces to be relined, on which the material is to adhere, evenly and thoroughly **with Softbase Primer**. In doing so, apply the primer approx. 1-2 mm beyond the border to provide security. Allow the primer to dry for approx. **1 minute**. Immediately close the vial with the primer after use to avoid evaporation. Saliva, moisture or monomer liquid prevent adhesion of Softbase to the denture.

#### 3. Preparing the cartridge for application

Place the cartridge in the dispenser 50 and lock it in place. Remove and discard the seal on the front end of the cartridge. Gently squeeze out the material until it flows evenly from both openings. Attach the mixing tip to the cartridge and squeeze out the material regularly.

**Notes:** The used mixing tip serves as a cap until the next application. Do not place the original cap back on the cartridge, as this makes it impossible to rule out cross contamination and thus unintentional partial curing.

### 4. Relining the denture

#### 4.1. Intraoral use

Apply Softbase evenly and extensively over the entire, previously prepared base surface of the denture. The applied quantity is based on the thickness indicated by the preparation and should be applied in excess. (Layer thickness at least 2 mm). It is important that all areas of the denture which come into contact with the gingiva are coated. Particular attention must be paid to avoiding air inclusions. The material application must be completed in **no more than 1:00 minute** (processing time at 23°C). (Including any smoothing performed using a clean spatula). Now immediately place the denture in the patient's mouth. The patient should close his/her mouth in normal firm occlusion. During the setting process, the patient performs a number of functional movements which ensure the subsequent hold of the denture. The setting process is complete after **approx. 4 minutes** and the denture can be taken out of the mouth for postprocessing as described in points 5 and 6.

#### 4.2. Extra oral use

Embed the denture to be relined in a conventional 2-part cuvette. In doing so, the row of teeth must be included such that the denture can be easily removed after opening. Take the denture out of the cuvette and then prepare as explained in points 1 and 2. Coat the plaster model using a conventional plaster-plastic insulation. Before applying the relining material, place the denture back in the cuvette. The prepared base surface points upward all around and is fully accessible. Using the dispenser, apply Softbase directly from the cartridge evenly and extensively on the denture area to be relined (layer thickness approx. 1 mm more than what was previously milled away). Close the cuvette with the counterweight and allow to stand with gentle pressure for **approx. 30 minutes** (at 23°C) under the press or **10 minutes** in the pressure pot at **40 to 45°C**. Then take the denture out of the cuvette and prepare and varnish it as described in points 5 and 6.

**Note:** Do not use with standard (condensation reaction type) silicone impression materials! Do not use any latex gloves since these can affect the setting reaction of the polyvinyl due to contamination.

### 5. Finishing the relined denture

Remove large amounts of excess soft liner using a sharp scalpel or scissors and then smooth the cut edges with a grinding tool suitable for silicones at 8000 rpm using gentle pressure.

### 6. Varnishing the denture

Perform final finishing of the post-processed areas by applying the **enclosed gloss varnish Softbase Varnish**. Prepare this gloss varnish, which is also based on addition-curing silicones, by thoroughly mixing 5 drops of each of the two components in the enclosed mixing palette, and apply with a brush. The varnished surface is set after a drying time of **5 minutes** at room temperature and the denture is ready for integration.

### 7. Cleaning recommendation for patients

The relined denture should be thoroughly cleaned several times per day using a toothbrush and neutral soap under running water; cleaning bath: max. 10 minutes.

### 8. Storage

Do not store at temperatures above 25°C. Do not use after the expiry date.

### Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico's sole obligation shall be repair or replacement of the Bisico product.

### Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

### Technical data

Mixing time:	N/A
Application time:	max. 1 min.
Retention time in the mouth:	4 min.
Total setting time:	5 min.
Polymerisation shrinkage:	< 0.1%
Compression set:	< 1.0%
Hardness (Shore A):	approx. 30

### For dental use only!

### Delivery form

**REF 73810**

4 ml-bottle Softbase Primer,  
50 ml cartridge Softbase, 6 mixing tips  
Softbase Varnish 5 ml bottle A, 5 ml bottle B  
1 brush holder, 10 disposable brushes, 1 mixing palette

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