Bisico® Glass Cem (Capsules)

Glass Ionomer Luting Cement

Bisico® Glass Cem is a newly developed glass ionomer luting cement that shows, even in warm and humid climate, an excellent shelf life.

Bisico® Glass Cem is resin free and shows very good biocompatibility.

Bisico® Glass Cem is designed for the following indications:

- cementing of crowns, bridges, inlays and onlays (metal, ceramic and zirconia)
- liner under composite fillings

Capsules are easily activated and the content of the capsule is easily ejected out with a special gun. Capsule mixing (mixing time 10 sec) is achieved by a high frequency mixer with about 4,300 oscillations/min such as Capmix. Application can be done directly from the Capsule.

1. Tooth Preparation

Prepare tooth in usual manner.

For optimal chemical bonding a **polyacrylic solution** can be applied for 20 sec to the bonding surfaces to remove the preparation smear layer. Rinse the **polyacrylic solution** thoroughly with water and dry gently. Do not desiccate. Pulp capping with **Bisico® Glass Cem** is contraindicated. To deep areas of possible pulpal exposure apply a small amount of calcium hydroxide liner. Prepare the carefully cleaned restoration according to the instructions of the manufacturer.

2. Activation and Mixing (see Instruction for Capsules)

Activate and mix the Capsules according to the information in the Capsule instruction. Mixing time for the Capsules is 10 seconds.

Attention: Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting which may impair or prevent application of the material. The material must be extruded within 10 seconds after the end of mix.

Bisico® Glass Cem sticks to metal instruments and should be washed off with cold water prior to setting.

3. Cementing

Apply the mixed material directly out of the Capsule onto the prepared bonding surface of the restoration and seat immediately (**working time 1:20 min** from start of mixing at 23°C).

Net setting time is about **4:30 min**.

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified.

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentin.

4. Conclusive Notes

The products are to be applied only by a dental professional in the manner as described in this instruction.

Do not use the products with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician.

Do not allow the liquid or the mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the mixture. In case of contact, immediately flush with water and seek medical treatment.

Storage

Store the products at room temperature. Do not use after expiry date.

Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico's sole obligation shall be repair or replacement of the Bisico product.

Limitation of Liability

Except where prohibited by law, Bisico, Bielefelder Dentalsilicone GmbH & Co.KG, will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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


Capsule before activation.


For activation of the Capsule press the plunger on a hard and plane surface to the end into the DirectCap.


Insert the Capsule into the DirectGun and **click once** to standardize.


Note: The plunger must be at the same level as the bottom of the capsule.

Insert the Capsule into a mixer (or an amalgamator), close lid and mix immediately (about 4300 oscillations / min). Mixing time: see corresponding recommendation for use

Remove the pin from the nozzle. If not, capsule can burst.

Insert the Capsule into the applicaton tool. Pull the lever 2 times (2 clicks) to prime the Capsule. Extrude the mixed material directly into the preparation. Unlock the gun (push button A) and remove the DirectCap.

Only with the application tool "DirectGun" the optimal amount of mixed material is guaranteed.