Algi-Press

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For dental professional use only.
 Federal (U.S.A.) Law restricts this device to sale or on the order of a dentist

 Bisico Bielefelder Dentsilicone GmbH & Co.KG
 Johanneswerkstraße 3
 33611 Bielefeld, Germany

Rev.: 65245.04-2024

Algi Press

Product description
 Algi-Press is a medium-viscosity, addition-curing silicone impression material, which combines the advantages of additioncuring silicones with those of alginates and can be used in place of typical alginate applications.

- Indications/Intended use**
 - Anatomical impressions
 - Opposing jaw impressions
 - Production of temporary restorations
 - Production of splints
 - Orthodontic applications

Patient target group
 Patients who require dental treatment or dentures

Users
 May only be used by qualified dental assistants and dentists.

Composition
 Polyvinyl polysiloxanes, fillers, pigments, additives, platinum compounds

Warnings
 We do not recommend using our silicones in combination with those from other manufacturers.
 Do not use impression materials with condensation-curing silicones, vinyl siloxane ethers, polyethers or polysulphide materials.
 Addition-curing silicones are chemically resistant and cannot be removed even by dry cleaning. The general hygiene measures must be observed.
 Only use Bisico silicones with the mixing tips provided by us (original Mixpac).

Contraindications
 If the patient is allergic or hypersensitive to any of the ingredients, this product should not be used or should only be used under the strict supervision of the treating dentist.

Hazard warnings
 Platinum compound, 4-methoxyphenol
 Normally, impression materials do not cause allergic reactions. However, an allergic reaction cannot be ruled out in sensitive patients.
 Do not use after the expiry date.

Interaction with other materials
 The use of latex gloves and various haemostatic agents can adversely affect the setting reaction of A-silicones. Once you have finished using these materials, these surfaces must be thoroughly cleaned to remove all residues.

Instructions for use
 In order to achieve optimum impressions, the product temperature should not deviate significantly from 23°C (74°F) so that the processing time and intraoral setting time are not affected.
 Only use components with the same batch designations.

Application
Preparation of the cartridge
 The use of Algi-Press requires a conventional mixing and dispensing machine. Please observe the instructions for use of these devices closely (on request, Bisico delivers the “Sympress” unit from Renfert, which is generally available on the market).
 Hold the cartridge with the sealed outlets facing upwards and remove the cap by grasping the tab and lifting upwards. Avoid exerting excessive force when doing so.
 Insert the cartridge into the output device as per the manufacturer's instructions and place a mixing tip over the outlets of the cartridge.
 Once the mixing tip is correctly positioned (the internal hexagon of the mixing tip must be correctly aligned with the hexagon of the drive shaft), push the yellow bayonet ring over the mixing tip as far as it will go and turn it clockwise (1/4 turn) to lock it securely in place.

Ensure that the drive hexagon engages in the mixing tip when the machine starts. The mixing spiral in the mixing tip should rotate when the machine is started. The tray can now be filled.

Notes
 When using a cartridge for the first time, discard the first 23 cm of material to adjust the levels of the two material components.
 Leave the used mixing tip on the cartridge as a cap until the next application. To remove the used mixing tip, release the bayonet ring by turning it anticlockwise and then pull of the mixing tip.

Bonding agent
 We recommend the use of Bisico Adhesive. Coat the impression tray with a thin layer, allow to dry for 1 minute and then cover the tray with a sufficient layer of Algi-Press.

Pouring the mouldings
 To obtain models of the highest surface quality, wait about 30 minutes before pouring the impression.
 The impression is permanently stable in storage and can be moulded several times.

Disinfection
 The impressions can be disinfected with all standard disinfectant solutions for impression silicones. The manufacturer's instructions should be followed in such cases.

Storage
 Bisico impression silicones should be stored at room temperature, approx. 25°C, in a dry place and protected from sunlight.
 Expiry date: see container

Disposal
 Impression silicones in cross-linked form may be disposed of with household waste in small quantities (e.g. sufficient for one impression). Impression material that has come into contact with the patient's mouth must be disposed of in accordance with local regulations.

Safety instructions
 In accordance with the EU Medical Devices Regulation, users are obliged to report serious incidents involving the medical devices to the manufacturer (stating the article number and batch) and the competent authority of the country in which the incident occurred.

Limitation of Liability
 In the event of material defects or obvious malfunctions, Bisico's liability shall be limited to the replacement of the affected material after inspection at our premises. Any further claims, in particular those resulting from the moulding methods used or from application errors, are excluded.

Warranty
 Bisico Bielefelder Dentsilicone GmbH & Co. KG warrants this product will be free from defects in material and manufacture. Bisico GmbH & Co. KG makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico GmbH's sole obligation shall be repair or replacement of the Bisico GmbH product.

For dental professional use only.
 Federal (U.S.A.) Law restricts this device to sale or on the order of a dentist

 Bisico Bielefelder Dentsilicone GmbH & Co.KG
 Johanneswerkstraße 3
 33611 Bielefeld, Germany

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