Nanofill flow

Intended purpose
Composite-based dental filling materials are used to build up and maintain tooth substance.

Product description
Nanofill flow is a light-curing, flowable, radiopaque composite with low viscosity.

Patient target group
All patients requiring tooth substance to be replaced or built up in order to restore or maintain dental function.

Users
Composite-based dental filling materials are used in dental practice or in a dental hospital by dental professionals.

Composition
Glass powder, diurethane dimethacrylate, Bis-GMA, tetramethylene methacrylate, splitter polymerisate
Filler content: 55% by weight (36% by volume) inorganic fillers (0.005 - 3 µm)

- Indications**
- Fissure sealing
 - Extended fissure sealing on molars and premolars
 - Fillings in Black's class V cavities (cervical caries, eroded areas in roots, wedge-shaped defects)
 - Minimally invasive fillings in Black's class I and II cavities in areas not exposed to severe occlusal loads
 - Minimally invasive fillings in Black's class III cavities
 - Restoring defects in enamel
 - Blocking out undercuts
 - Minimal adjustments to the contours and shade of the enamel

Contraindications
If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist. Linings containing eugenol are contraindicated.

Warnings
Contains: tetramethylene dimethacrylate, di-urethane dimethacrylate
Warning: May cause allergic skin reactions. Harmful to aquatic life with long lasting effects.

Safety instructions
Avoid breathing vapours/spray. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/medical attention.

Interactions with other materials
Phenolic substances (such as eugenol) inhibit polymerisation. Therefore, do not use any lining materials (such as zinc oxide-eugenol cements) which contain substances of this type. The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

Application
Pre-treatment
Prior to the treatment, clean the hard tooth tissue with a fluoride-free polishing paste. Select colours when still wet using the Vita® Shade Guide.

1. Cavity preparation
Prepare the cavity minimally invasively as generally required for adhesive techniques. All enamel margins in the anterior region must be beveled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.

2. Pulp protection / Cavity liner
If an enamel-dentine adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.

3. Approximal contact areas
For cavities with approximal sections, place a transparent matrix and fix it in place.

4. Adhesive system
Etch and bond according to the manufacturer's instructions. We recommend the use of E-Bond LC.

5. Applying Nanofill flow
Introduce Nanofill flow in thin layers (max. 2 mm) directly into the cavity. Use the curved delivery tips included with the kit for application. For hygiene reasons, tips are intended only for single use. Avoid introducing air bubbles during application with the tip. Ensure the prepared hard tooth tissue is thoroughly wetted.

6. Curing
Cure each layer after coating for 40 seconds using a commercially available polymerisation unit. The light guide must be held as close as possible to the surface of the filling. A dispersion layer forms on the surface during polymerisation which should not be touched or removed. It is used to bond the following composite layers.

7. Finishing
Nanofill flow can be finished and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishers, and polishing brushes can be used for finishing.

Please note
• When placing time consuming restorations, to prevent the composite curing prematurely the dental light should be moved away from the site temporarily or the composite covered with foil impervious to light.

- For hygienic reasons, the curved application tips supplied with the material must only be used once!
- Use a light polymerization system with an emission range of 350 - 500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions.

Light intensity for polymerization	≥ 650 mW/cm²
Wave length for polymerization	350 - 500 nm
Polymerization time	40 sec.
Maximum layer thickness	2 mm

Troubleshooting	Trouble	Cause	Remedy
Composite does not cure		Luminous intensity of the polymerization unit insufficient	Check luminous intensity; replace light source, if necessary
		Emitted spectral range of the polymerization unit insufficient	Consult manufacturer of polymerization unit; recommended spectral range: 350 - 500 nm
Composite seems to be too hard/ firm inside the syringe		Material was stored at temperatures below 10°C for a longer period of time	Let composite reach room temperature before use
		Syringe was not closed tightly which caused part of the material to cure	Close syringe correctly with the cap after each use
Composite does not cure sufficiently		Layer thickness per polymerization cycle too high	Keep to max. layer thickness of 2 mm
Restoration seems too yellow when compared to color reference		Insufficient polymerization of the composite layers	Repeat polymerization cycle several times, for a minimum of 40 sec.

Information on storage and handling
Store at 10 - 25°C (50 - 77°F). Close syringe tightly immediately after use. The material should be at room temperature before use. Retract the syringe spindle slightly after use to prevent the discharge opening from becoming clogged.

Shelf life
The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

Side effects
With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

Disposal
Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

*Vita is a registered trade mark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.

Delivery Form		
Nanofill flow, A1		74062
Syringe with 2 g, 5 Appli-Tips		
Nanofill flow, A2		74061
Syringe with 2 g, 5 Appli-Tips		
Nanofill flow, A3		74060
Syringe with 2 g, 5 Appli-Tips		
Nanofill flow, A3.5		74070
Syringe with 2 g, 5 Appli-Tips		

Bielefelder Dentsilsilicone GmbH & Co. KG
Johanneswerkstraße 3
33611 Bielefeld
Tel.: + 49 (0) 521 8 01 68 00
Fax: + 49 (0) 521 8 01 68 01

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*Vita ist eine eingetragene Marke der Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.
MANI MEDICAL GERMANY GmbH