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Instructions for use **EN**

Provi Temp K
Temporary Crown & Bridge Material

Provi Temp K is an automatically mixed two component material in cartridges based on multifunctional methacrylic esters. **Provi Temp K** is a very easy and comfortable material for creating short- and long-term temporary crowns or bridges, Inlays, onlays and veneers. Because of its flexibility, the material is especially suitable for longer bridge spans. **Provi Temp K** shows fluorescence in UV-light.

Provi Temp K is free of methyl methacrylate. Its temperature derived while curing is lower than 40 ° C. As a temporary crown or bridge it protects the prepared teeth against external influences and preserves the occlusion. It shows increasing transverse strength, good abrasion resistance and low polymerization shrinkage. It fits perfect. Showing good polishability, good color stability and fluorescence it looks perfect.

And if **Provi Temp K** breaks - which may happen with acrylics - it is easy to repair. Cured material can be repaired with any natural or light cure composite building up again a chemical link. In nearly every property **Provi Temp K** is superior to the common materials based on monomer/polymer systems. It allows thin edges and production of longer bridge spans.

Application

1. Impression taking

Before preparing of stumps or before a planned extraction of a tooth a situation impression is made by addition curing silicones (long storage stability) or by alginates. Interdental areas are carved out, in molar areas with teeth absent it may be necessary to cut a groove in the impression between the abutments to create a bridge-like connection between the tooth units.

Note: In the silicone impression block out undercuts and if necessary cut grooves into the impression.

2. Preparing of the cartridge

The turn cap of the cartridge is taken o and thrown away (do not use it again!). It is substituted by the supplied mixing cannula by attaching it with a sideways turning movement, locking the cannula in place. The application gun is loaded with the prepared cartridge and is ready for application.

Note: The initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix. Leave the used mixing cannula on the cartridge. It serves as a cap.

3. Application

Provi Temp K is automatically mixed when dispensed with slight and even pressure directly into the situation impression made before. Filling should occur from bottom upward to prevent voids.

4. Forming of the temporary crowns or bridges

- Load the situation impression with Provi Temp K.
- Seat the impression onto the prepared areas of the teeth. Step a and b must be done within 50 seconds (working time at 23°C).
- after 1-2 minutes (setting time in mouth) the material shows a hardened but still elastic condition and can be removed from the teeth (**alternatively: 3-4 minutes setting time on the model**)

Note: The setting reaction has to be checked with excess material intra-orally (e.g. with a dental probe) resp. on the model. The oral temperature has a significant eect on the setting reaction and the temporary can only be removed without destruction during the elastic state.

5. Post curing and finishing

After removal of the temporary crown or bridge from the impression (or from the tooth preparation) excess material and proximal undercuts are removed. Preferably, the temporary is then post cured in warm (45° - 55° C) water (e.g. in a hot cure polymerization device). If this is not possible or not desired, after 6 minutes from beginning of application the temporary acrylic can be worked out with rotary instruments and can be polished to high gloss. Do not breathe polishing dust; use suitable mouth protective device or aspiration!

Note: The oxygen inhibited smear layer on the surface caused by the acrylic system can not be fully avoided and should be removed before working out. It can easily be removed by alcohol or other suitable solvents.

6. Cementing of the temporary

Provi Temp K temporary crowns or bridges preferably should be cemented with an eugenol-free temporary cement (e.g. Provi Bond K). In case of using eugenol-containing cements it should be considered that later used acrylate based cements could be hindered in hardening (inhibition of curing reaction by traces of remained eugenol).

7. Repairs

Provi Temp K temporaries show high mechanical strength. However, if a temporary breaks the following procedures are recommended:

a) Fracture shortly after production

Both ends of the fracture are refixed with freshly extruded Provi Temp K. Fine flowing composite can be applied alternatively.

b) Fracture of longer existing temporary:

The areas of fracture are cleaned and roughened and provided with some mechanical retentions. Thus prepared, the fractured areas can be joint with freshly extruded Provi Temp K to enhance the curing place for some minutes in warm water. A flowable light cure composite (e.g. Bisico Nano fill flow) can be applied alternatively. To ensure a complete cure, do not apply Nano fill flow in layers exceeding 2 mm. Light cure each layer for 40 seconds..

8. Special hints

Non hardened Provi Temp K can easily be removed with alcohol or other suitable solvents.

9. Precautions

Provi Temp K is free of methyl methacrylate but contains other acrylates. With susceptible patients, sensitization to Provi Temp K cannot be excluded. Provi Temp K should not be used any more, if allergic reactions are observed. Do not use for patients with allergic reactions against acrylates. Avoid contact with skin, mucous membrane and eyes If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.

10. Storage


Do not store above 25 ° C! Provi Temp K shall not be used after expiry date!

Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico's sole obligation shall be repair or replacement of the Bisico product.

Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Time Table of application and curing of Provi Temp K				
0:00	0:50 min.	approx. 3 min	approx. 4:30 min.	
				
temperature of the material [°C]	approx. 23	approx. 37	approx. 50	approx. 23
	If the completion of setting takes place at room temperature, trimming / polishing could be done not before after 6 minutes (from start of mixing)			

Keep away from children! For dental use only!	
Delivery form Provi Temp K, A2 1 cartridge of 50 ml, 15 mixing tips light blue 4:1	REF 73540
Provi Temp K, A2 1 mini syringe of 5 ml, 5 mixing tips brown 4:1	REF 73541
Provi Temp K, A1 1 cartridge of 50 ml, 15 mixing tips light blue 4:1	REF 73560
Provi Temp K, A3 1 cartridge of 50 ml, 15 mixing tips light blue 4:1	REF 73570
Provi Temp K, bleach 1 cartridge of 50 ml, 15 mixing tips light blue 4:1	REF 73590

  0482	
  0482	
  0482	
  0482	
  0482	

  0482	
  0482	
  0482	
  0482	
  0482	

  0482	
  0482	
  0482	
  0482	
  0482	

  0482	
  0482	
  0482	
  0482	
  0482	

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Bisco GmbH