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Micro Esthetic flow-viscous/ Micro Esthetic Gingiva

Produktbeschreibung

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Micro Esthetic flow-viscous/ Micro Esthetic Gingiva

Product Description

Micro Esthetic flow-viscous/ Micro Esthetic Gingiva is a light-curing, flowable, highly radiopaque (210% Ad) composite with high viscosity. The guidelines and specifications of DIN EN ISO 4049 apply.

Micro Esthetic flow-viscous/ Micro Esthetic Gingiva composite is available in syringes.

Intended purpose

Composite-based dental filling materials are used to build up and maintain tooth substance.

Patient target group

All patients requiring tooth substance to be replaced or built up in order to restore or maintain dental function.

Users

For use in the dental practice by dental professionals

Composition

Glass powder, diurethane dimethacrylate, silicon dioxide, tetramethylene dimethacrylate

Filler content: 77% by weight (57% by volume) inorganic fillers (0.005 - 40 µm)

Indications

- Fissure sealing
- Extended fissure sealing on molars and premolars
- Fillings of class V according to Black (cervical caries, root erosions, wedge-shaped defects)
- Minimally invasive fillings of classes I, II and III according to Black
- Correction of dental enamel defects
- Blocking out of undercuts
- Minor shape and colour corrections on the enamel

Contraindications

If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/ dentist.

Hazard and safety information

Contains: tetramethylene dimethacrylate, di-urethane dimethacrylate.

Warning: May cause an allergic skin reaction. Avoid breathing vapours/spray. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/medical attention. Avoid release to the Environment.

Interactions with other materials

Phenolic substances (such as eugenol) inhibit polymerisation. Therefore, do not use any lining materials (such as zinc oxide-eugenol cements) which contain substances of this type.

The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

1. Application

Prior to the treatment, clean the hard tooth tissue with a fluoride-free polishing paste. Select colours when still wet using the Vita** Shade Guide.

1.2 Cavity preparation

Preparation of the cavity with protection of the hard tooth tissue according to the general rules of adhesive technique. In the region of the anterior teeth, all enamel edges must be bevelled. In the region of the posterior teeth, by contrast, do not bevel the edges and avoid feather edges. Then clean the cavity with water spray, remove all residues, and dry it. Isolation is necessary. The use of a rubber dam is recommended.

In case of non carious defects the tooth surface has to be roughening before using Micro Esthetic Gingiva.

1.3 Pulp protection/lining

Lining is not necessary if an enamel-dentin adhesive is used. If there are very deep cavities near the pulp, cover corresponding areas with a calcium hydroxide (e.g. BisiCAL) preparation.

1.4 Approximal contact areas

For cavities with approximal sections, a transparent matrix must be created and fixed in place.

1.5 Adhesive system

Etch (z. B. Best Etch) and bond (z. B. Bond LC) according to the manufacturer's instructions. For Micro Esthetic Gingiva we recommend a selfetching adhesive, e.g. E-Bond LC. Please use the adhesive according to the manufacturer's instructions.

1.6 Application

Rotate the delivery tip clockwise until it locks in place on the syringe and dispense Micro Esthetic flow-viscous/ Micro Esthetic Gingiva in thin layers (max. 2 mm) directly into the cavity. Use the curved delivery tips included with the kit for application. For hygiene reasons, tips are intended only for single use. Avoid introducing air bubbles during application with the tip. Ensure the prepared hard tooth tissue is thoroughly wetted.

1.7 Curing

Cure each layer after coating for 40 seconds using a commercially available polymerisation unit. The light guide must be held as close as possible to the surface of the filling. A dispersion layer forms on the surface during polymerisation which should not be touched or removed. It is used to bond the following composite layers.

1.8 Finishing

Micro Esthetic flow-viscous/ Micro Esthetic Gingiva can be finished and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishers, and polishing brushes can be used for finishing. Check the occlusion and articulation and grind so that no early contacts or undesirable paths of articulation remain on the surface of the restoration.

Special information

- In case of time-consuming restorations, the surgical light should be temporarily moved further away from the work area to prevent premature curing of the composite or the material should be covered with an opaque foil.
- The curved delivery tips for the application which are supplied with the kit are intended for single use only since contamination of the materials and spread of germs cannot be ruled out.
- A light polymerisation unit with an emission spectrum in the range of 350 - 500 nm is to be used for polymerisation. The physical properties required are achieved only with properly working lamps. For this reason, it is necessary to regularly check the light intensity according to the manufacturer's instructions.

Light intensity for the curing	≥ 650 mW / cm ²
Wavelength for the curing	350 - 500 nm
Curing time	40 sec.

Troubleshooting

Fehler	Ursache	Abhilfe
Composite does not cure	Light output of the polymerisation lamp is inadequate	Check the light output Clean light guide if soiled. Replace the light source, if necessary
	Emitted wavelength range of the polymerisation lamp is inadequate	Consult manufacturer of the polymerisation lamp. Recommended wavelength range: 350 - 500 nm
Composite appears to be too hard and solid in the syringe	Material stored for a longer period of time at temperatures <10°C	Allow composite to warm up to room temperature before use
	Syringe not closed correctly, composite polymerised	Close the syringe correctly with the cap each time after removing composite
Composite does not correctly cure (dark or opaque shades)	Composite layer too thick per curing cycle	Maintain max. layer thickness of 2.0 mm per layer
Restoration appears too yellow in comparison to the colour reference	Inadequate polymerisation of the composite layering	Repeat exposure cycle several times, min. 40 sec.

Information on storage and handling

Store at 10 - 25°C (50 - 77°F). Close syringes tightly immediately after use. The material should be at room temperature before use. Retract the syringe spindle slightly after use to prevent the discharge opening from becoming clogged.

Shelf life

The maximum shelf life is printed on the label of each syringe or directly on the complete. Do not use after the expiry date.

Side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

Disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

*Vita is a registered trade mark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad

Bisico Bielefelder Dental silicone GmbH & Co.KG
Johanneswerkstraße 3
33611 Bielefeld, Germany

MANI MEDICAL GERMANY GmbH
Hertha-Sponer-Straße 2
61191 Rosbach v.d. Höhe, Germany

Rev.: 700170 v. 20230901

Rev.: 700170 v. 20230901