Nanofill unique flow

Intended purpose

Composite-based dental filling materials are used to build up and maintain tooth substance.

Product description

Nanofill unique flow is a light-curing, flowable, highly radiopaque (210% AI) composite with high viscosity. The guidelines and specifications of DIN EN ISO 4049 apply.

Nanofill unique flow is available in compules. The compules are for single use. Please do not reuse them since contamination and germ formation otherwise cannot be excluded.

Patient target group

All patients requiring tooth substance to be replaced or built up in order to restore or maintain dental function.

Users

For use in the dental practice by dental professionals.

Composition

Glass powder, diurethane dimethacrylate, silicon dioxide, tetramethylene dimethacrylate

Filler content: 77% by weight (57% by volume) inorganic fillers (0.005 - 40 µm)

Indications

- Fissure sealing
- Extended fissure sealing on molars and premolars
- Fillings in Black's class V cavities (cervical caries, eroded areas in roots, wedge-shaped defects)
- Minimally invasive fillings in Black's class I, II, III
- Correction of dental enamel defects
- Blocking out undercuts
- Minor shape corrections on the enamel

Contraindications

If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist.

Hazard and safety information

Contains: tetramethylene dimethacrylate, di-urethane dimethacrylate.

Warning: May cause an allergic skin reaction. Avoid breathing vapours/spray. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/medical attention. Avoid release to the Environment.

Interactions with other materials

Phenolic substances (such as eugenol) inhibit polymerisation. Therefore, do not use any lining materials (such as zinc oxide-eugenol cements) which contain substances of this type. The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

Application

Prior to the treatment, clean the hard tooth tissue with a fluoride-free polishing paste.

1. Cavity preparation

Preparation of the cavity with protection of the hard tooth tissue according to the general rules of adhesive technique. In the region of the anterior teeth, all enamel edges must be bevelled. In the region of the posterior teeth, by contrast, do not bevel the edges and avoid feather edges. Then clean the cavity with water spray, remove all residues, and dry it. Isolation is necessary. The use of a rubber dam is recommended.

2. Pulp protection / Cavity liner

Lining is not necessary if an enamel-dentin adhesive is used. If there are very deep cavities near the pulp, cover corresponding areas with a calcium hydroxide preparation.

3. Approximal contact areas

For cavities with approximal sections, place a transparent matrix and fix it in place.

4. Adhesive system

Etch and bond according to the manufacturer's instructions. We recommend the use of E-Bond LC.

5. Applying Nanofill unique flow compules

Insert the compule in the dispenser. Remove the cap. Secure the compule such that the opening is aligned at the correct angle for application into the cavity. Introduce the material into the cavity using slow, even pressure. Do not use excessive force! The layer thickness must not exceed 2 mm. When finished, pull back the spindle to remove the compule from the dispenser. Then remove the compule.

Note: For hygiene reasons, compules are intended only for single use.

6. Curing

Cure each layer after coating for 40 seconds using a commercially available polymerisation unit. The light guide must be held as close as possible to the surface of the filling. A dispersion layer forms on the surface during polymerisation which should not be touched or removed. It is used to bond the following composite layers.

7. Finishing

Nanofill flow can be finished and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishers, and polishing brushes can be used for finishing. Check the occlusion and articulation and grind so that no early contacts or undesirable paths of articulation remain on the surface of the restoration.

Please note

• In case of time-consuming restorations, the surgical light should be temporarily moved further away from the work area to prevent premature curing of the composite or the material should be covered with an opaque foil.

• A light polymerisation unit with an emission spectrum in the range of 350 - 500 nm is to be used for polymerisation. The physical properties required are achieved only with properly working lamps. For this reason, it is necessary to regularly check the light intensity according to the manufacturer's instructions.

Light intensity for polymerization	≥ 650 mW/cm ²
Wave length for polymerization	350 - 500 nm
Polymerization time	40 sec.

Troubleshooting

Trouble	Cause	Remedy
Composite does not cure	Light output of the polymerisation lamp is inadequate	Check luminous intensity; replace light source, if necessary
	Emitted wavelength range of the polymerisation lamp is inadequate	Consult manufacturer of polymerization unit; recommended spectral range: 350 - 500 nm
Composite does not cure sufficiently	Composite layer too thick per curing cycle	Keep to max. layer thickness of 2 mm
Restoration seems too yellow when compared to color reference	Inadequate polymerisation of the composite layering	Repeat polymerization cycle several times, for a minimum of 40 sec.

Information on storage and handling

Store at 10 - 25°C (50 - 77°F). The material should be at room temperature before use.

Shelf life

The maximum shelf life is printed on each compule. Do not use after the expiry date.

Side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

Disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

Bielefelder Dental silicone GmbH & Co. KG
Johanneswerkstraße 3
33611 Bielefeld
Tel: + 49 (0) 521 8 01 68 00
Fax: + 49 (0) 521 8 01 68 01

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MANI MEDICAL GERMANY GmbH