Produktbeschreibung

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BisiMTA
Endodontic Repair Cement

Product description

BisiMTA is a biocompatible endodontic repair cement. **BisiMTA** powder is consisting of very fine hydrophilic particles of several mineral oxides. After contact with **BisiMTA** liquid it forms a gel that hardens to an impermeable barrier.

Indications/Intended use

- Repair of root perforations (during root canal therapy)
- Root end filling (retrograde)
- Pulp capping
- Apexification / Root end filling (orthograde)
- Repair of internal resorption
- Pulpotomy

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

None known.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Side effects

None known

Application

1. Dispensing and Mixing

The powder/liquid ratio is 2.6/1.0. This can be obtained by mixing 1 level (blue) scoop of powder and 2 drops of liquid.

If a thinner or firmer consistency is desired, the mixing ration can be modified slightly:

Mixing ratio (by weight)	2:1	2,6:1	3:1
VWorking time (at 23 °C/74 °F)	3:00 min	2:00 min	1:00 min

For root end filling (orthograde) (see 5.4) the recommended mixing ratio is 2:1 (by weight). This can be obtained by mixing 3 level (blue) scoops of powder and 8 drops of liquid. The mixed material gives a sufficient amount for apexification that can be applied optimally with a suitable application device into the root canal.

For mixing of **BisiMTA** use a mixing pad that is impervious to water or a glass block of suitable dimension.

For accurate dispensing of **BisiMTA** powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder with the mixing spatula and carry it onto the mixing pad.

For dispensing of **BisiMTA** liquid turn the bottle vertically with the tip about 5cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop at a time. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discount under-sized drops that contain bubbles and are obviously not full-sized. Discount over-sized drops, usually resulting from holding the bottle too close to the mixing pad or squeezing the bottle too hard and/or for too long.

Use a small spatula to rapidly mix all the cement powder in portions into the liquid. The mixed cement should be thixotropic and have a homogeneous consistency. Total mixing time is 30 seconds.

If desired, a more rigid consistency can be achieved by adding some more powder to the mixture, a more creamy consistency is attained by adding some liquid.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

2.Application

2.1. Repair of root perforations

Place rubber dam and clean the root canal system using intracanal instruments and NaOCl solution. Dry the root canal with paper points and isolate the perforation.

Fill the root canal apical of the perforation completely with a suitable root canal filling material.

Mix **BisiMTA** as described under point 1.

Apply **BisiMTA** with suitable instruments into the perforation site and condense it.

Check the position of **BisiMTA** in the root canal by an X-ray. If an adequate barrier has not been created, rinse **BisiMTA** out of the canal and repeat the procedure.

Remove excess moisture with a cotton pellet or a paper point.

Place a cotton pellet in the access to the root canal and apply a temporary filling material. Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling.

Both options can be done at the earliest 5 minutes after placement of the **BisiMTA**. **BisiMTA** repair material remains as a permanent part of the root canal filling.

2.2. Root End Filling (retrograde)

Create a surgical access to the root-end and resect the root. Prepare an apical cavity to a depth of 3-5 mm.

Isolate the area and dry the root end cavity with paper points. Achieve hemostasis with suitable methods.

Mix **BisiMTA** as described under point 1.

Apply **BisiMTA** with suitable instruments and condense it using a small plugger. Remove excess cement and clean the surface of the root with a moist piece of gauze. Confirm placement of the **BisiMTA** repair material with an X-ray. **BisiMTA** remains as a permanent part of the root canal filling.

2.3. Pulp Capping

Place rubber dam and prepare the cavity. Rinse the cavity and exposed pulpal areas with a suitable disinfectant.

Mix **BisiMTA** as described under point 1.

With a suitable instrument apply a small amount of **BisiMTA** over the exposed pulp and remove excess moisture with a cotton pellet.

At the earliest 5 minutes after application of **BisiMTA** place a small amount of a flow-able light cure liner and light cure.

Etch the remaining cavity walls according to the total-etch-technique with ETCHING GEL (e.g. Bisico Best Etch) and apply a suitable bonding agent (e.g. Bond LC) according to the corresponding instructions.

Place a light cure composite (e.g. Nanofill n) according to the instructions and light cure.

Check pulp vitality and status regularly.

2.4. Apexification / Root End Filling (orthograde)

Place rubber dam and clean the root canal system using intracanal instruments and irrigate with NaOCl. Dry the root canal with paper points.

For disinfection place calcium hydroxide paste in the root canal and seal the access opening with a temporary filling material.

After one week remove the calcium hydroxide paste from the root canal.

Mix **BisiMTA** as described under point 1.

With a suitable instrument apply a small amount of **BisiMTA** into the apical region and condense it. Create a 3 - 5 mm barrier of **BisiMTA**.

Check the position of **BisiMTA** by an X-ray. If an adequate barrier has not been created, rinse **BisiMTA** out of the canal and repeat the procedure.

At the earliest 5 minutes after application of **BisiMTA** place a cotton pellet in the access to the root canal and apply a temporary filling material. Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling.

Additional Notes/Warnings

- In the first hour after application handle the placed MTA cement carefully.
- **BisiMTA** can cause discoloration.
- Keep away from children

Composition

MTA, bismuth oxide

Storage

Store **BisiMTA** protected from moisture at 10-25 °C (50-77 °F). Do not use after expiry date.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED - <https://ec.europa.eu/tools/eudamed>).

Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico GmbH's sole obligation shall be repair or replacement of the product.

Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.