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Nanofill - lichthärtendes Nano-Komposit

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Bisico Nanofill

Instructions for use **EN**

Nanofill is a light-curing nano-composite for the adhesive filling technique. It contains an ultrafine, radiopaque glass filler. The composition with ultrafine filler leads to exceptionally homogenous restorations. An optimal adaptation of the color is the result of an optimally adjusted chameleon effect. The guidelines and requirements of DIN EN ISO 4049 apply. Nanofill should be used in connection with the system components Best Etch and Bond LC or E-Bond LC. Nanofill is available in syringes and compules. The compules are for single use only. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation.

Composition

Glass powder, silicon dioxide, diurethane dimethacrylate, Bis-GMA, tetramethylene dimethacrylate, prepolymer

Filler content 83.5 % by weight (66.5 % by volume)
Anorganic filler (28 nm- 40 µm)

Indications

- Direct anterior and posterior tooth restorations in Black's classes I, II, III, IV and V.
- Indirect restorations such as inlays, onlays and veneers
- Extended fissure sealing on molars and premolars
- Building up stumps
- Splinting of loosened teeth
- Corrections of shape and color to enhance aesthetics

Available colors

4 g syringe: A1, A2, A3, A3.5, A4, B1, B2, B3, C (Clear), EW, OA, bleach

Application

Preparation of the natural tooth

Before starting the treatment, clean the hard tissue with fluoride-free polishing paste. Choose the color with help of the Vita*-shade guide while the tooth is still moist.

Cavity preparation

Gently prepare the hard tissue according to the rules of the adhesive technique. When working on anterior teeth, bevel all enamel edges. Do not bevel the edges when working on posterior teeth. Avoid feather-edged margins. Next, clean the cavity with water spray, remove all residue and dry. Complete drying is necessary. We recommend the use of a kofferdam.

Pulp protection / liner

When using an enamel-dentin adhesive the liner can be foregone. In case of very deep cavities near the pulp, line with a calcium hydroxide compound accordingly.

Design of approximal contacts

In cavities with approximal sections, set in a transparent matrix and fix.

Production of inlays, onlays and veneers

At the lab, pour an extra-hard plaster into the impression. Once the die has hardened, remove the impression from the die. Block out undercuts and insulate the die with an oil-free insulator. Build up the inlay on the die layer by layer. First, build up approximal and deep occlusal areas. Each layer should have a maximum thickness of 2 mm.

Polymerize with a commercial polymerization system. Lift off the completed restoration from the die, finish and polish to a high gloss. Clean the restoration thoroughly with water and soap. Rinse with air/ water spray and dry.

Inserting inlays, onlays and veneers

Remove the temporary restoration and clean the cavity. Set in a kofferdam, then clean and dry the surface of the tooth. Check the restoration's fit with slight pressure. Avoid rough handling. If necessary, improve the fit by grinding the inside surface. Do not check the occlusion in this first fitting, as the restoration could break.

Etch (e. g. Best Etch) and bond (e. g. Bond LC or E-Bond LC) according to manufacturer's instructions.

Fixing the restoration

The restoration is fixed with a commercially available, dual-curing fixing composite. Please adhere to the manufacturer's instructions.

Special instructions

- The working range under a surgical light is 2 min
- In case of extensive restorations, the surgical light should be moved away from the working area temporarily to avoid premature curing of the composite. Alternatively, the material can be covered with a light-tight foil.
- Use a light polymerization system with an emission range of 350 - 500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions.

Light intensity for curing	≥ 1200 mW/cm ²
Wavelength for curing	350 – 500 nm
Curing time	20 sec.

Contains teramethylene

Warning: May cause an allergic skin reaction. Avoid breathing vapours/ spray.

Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. In order to avoid a possible pulp reaction, always prepare a liner in cases with exposed dentin (e. g. compound containing calcium hydroxide).

Contraindications / interactions:

If a patient has known allergies against or hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product. Phenolic substances (e. g. eugenol) inhibit polymerization. Therefore, these materials (e. g. zinc oxide eugenol cements) must not be used as liners.

Storage

Store at temperatures between 10 °C and 25 °C (50 °F to 77 °F) . Avoid direct sunlight. Screw the cap back onto the syringe tightly after each use. Let the material reach room temperature before use. Withdraw the plunger slightly after use to keep the outlet from becoming plugged. Do not use after the expiration date (see label on the syringe). Only for use in dentistry. Keep out of children's reach. This product has been developed for the specific use illustrated above. Only process as described in these instructions. The manufacturer will not be held liable for any damages that result from improper use or improper processing. We recommend to forego the use of syringe heaters.

Adhesive system

Etch (e. g. Best Etch) and bond (e. g. Bond LC or E-Bond LC) according to manufacturer's instructions.

Application of the composite

Take the required amount of composite out of the application tip. Fill the cavity with the material and shape as needed, using conventional metal instruments.

A layer thickness of 2 mm must not be exceeded. Due to the oxygen in the ambient air, a thin dispersion layer will remain on top of each layer. This dispersion layer forms the chemical bond between the layers and must not be touched or contaminated with moisture.

Polymerization

The polymerization time is 20 sec. per layer for all colors in a halogen polymer- ization system or an LED polymerization lamp or 2 x 3 sec. in a plasma poly- merization system. Hold the light-guide as close to the surface of the filling as possible. Fillings with several surfaces should be polymerized from the direction of each surface. Due to the oxygen in the ambient air, a thin dispersion layer will remain on top of each layer. This dispersion layer forms the chemical bond between the layers and must not be touched or contaminated with moisture.

Finishing

Nanofill can be finished and polished immediately after polymerization. Suitable are finishing diamonds, flexible separating discs, silicone polishers and polishing brushes. Check occlusion and functional movements, correct if necessary. Finally, polish with suitable polishing pastes.

Indirect method

Cavity preparation:

A tooth-conserving preparation with only little divergence in the cavity walls should be preferred. All inner edges and angles must be rounded, avoid feathery margins. Design a planar cervical shoulder, do not bevel it. Block out any unavoidable undercuts with glass ionomer cement. For the preparation, use slightly cone-shaped diamond grinders with rounded edges. Cover dentin close to the pulp with a thin layer of a calcium hydroxide compound (e.g. BisiCAL). Any liner materials that contain eugenol are contra-indicated.

Impression and temporary restoration

After taking the impression, a temporary resin restoration is produced. Fix temporaries only with eugenol-free cement.

* Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.

Trouble shooting

Trouble	Cause	Remedy
Composite does not cure.	Light output of the polymerization lamp is insufficient.	Check the light output and, if necessary, replace the light source.
	Emitted wavelength range of the polymerization lamp is insufficient.	Consult the manufacturer of the polymerization lamp. Recommend wavelength range 350 - 500 nm.
Composite inside the syringe is sticky and soft; clear liquid separates from the material inside the syringe.	Material was stored at temperatures ≥ 25 °C (≥ 77 °F) for a longer period of time.	Please note the storage temperatures: Store at temperatures between 10 °C (50 °F) and 25 °C (77 °F). Short-duration storage in a refrigerator.
	Material was stored in the syringe heater for too long.	Never leave a syringe longer than one hour per use in a syringe heater.
Composite appears to be too compact and hard inside the syringe.	Material was stored at temperatures < 10 °C (< 50 °F) for a longer period of time.	Let composite reach room temperature before use; optionally use a syringe heater for a brief time.
	Syringe was not closed correctly, composite has started to polymerize.	Each time after removing material from the syringe, correctly cap and close the syringe.
Inlay/onlay does not stay in place after insertion.	Restoration appears too opaque to be fixed with a purely light-polymerizing composite.	Use dual-polymerizing fixing composite.
Composite does not cure all the way through (dark or opaque colors).	Layers per polymerization cycle were too thick.	Do not exceed max. layer thickness of 2.0 mm per layer.
Restoration has a yellowish tint when compared to the color reference.	Insufficient polymerization of the composite layers.	Repeat polymerization cycle several times; at least 20 sec.

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Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, how ever, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. In order to avoid a possible pulp reaction, always prepare a liner in cases with exposed dentin (e. g. compound containing calcium hydroxide).

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* Vita è un marchio registrato di Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.