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GLASS CEM - HandMix

Glasionomer Befestigungszement

Produktbeschreibung

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Glass Ionomer Luting Cement

Product description

GLASS CEM is a new developed glass ionomer restorative cement that shows, even in warm and humid climate, an excellent shelf life.

GLASS CEM is resin-free and shows very good biocompatibility.

Indications/intended use

- cementing of crowns, bridges, inlays and onlays (metal, ceramic and zirconia)
- liner under composite fillings

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

- Pulp capping
- Allergic reactions to glass ionomer cements

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Application

1. Tooth Preparation

Prepare tooth in usual manner.

For optimal chemical bonding a dentine conditioner can be applied according manufacturer's instructions for use.

Pulp capping with GLASS CEM is contraindicated. To deep areas of possible pulpal exposure apply a small amount of BisiCAL.

Prepare the carefully cleaned restoration according to the instructions of the manufacturer.

2. Dispensing and Mixing

The powder/liquid ratio to achieve a suitable consistency is 1.2 / 1.0. This can be obtained by mixing 1 level scoop of powder and 2 drops of liquid.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

For accurate dispensing of GLASS CEM Powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder by using the collar at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon with the inside wall of the bottle.

For dispensing of GLASS CEM Liquid turn the bottle vertically with the tip about 5cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense the drops. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discard drops that contain bubbles and are obviously not full-sized.

Use a small spatula to rapidly mix half the powder into the liquid.

Then mix the remaining powder into the mixed paste to achieve a glossy appearance. Total mixing time is 30 sec.

After use, close both powder and liquid bottles.

3. Cementing

Mix the required amount of GLASS CEM and apply it to the prepared bonding surface of the restoration and seat immediately (working time 1:30 min from start of mixing at 23°C).

Net-setting time is about 4:30 min.

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified (approx. 5:00 min).

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentin.

Storage

Store the products at room temperature (4 - 25 °C / 39 - 77 °F). Do not use after expiry date.

Additional Notes/Warnings

- Do not use the products with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician.
- Do not allow the liquid or the mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.
- Avoid eye contact of the mixture. In case of contact, immediately flush with water and seek medical treatment.
- Do not mix the powder or liquid of GLASS CEM with any other restorative cement components.
- Keep away from children!

Composition

Dental glass, polyalkenoate acids.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED - <https://ec.europa.eu/tools/eudamed>).

Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Bisico GmbH's sole obligation shall be repair or replacement of the Bisico GmbH product.

Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.