Micro Esthetic

bisico®

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Light Cure Microhybrid Composite

bisico® Instructions for Use EN

Product description
Micro Esthetic is a visible light cure radioaque microhybrid composite for anterior and posterior restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0.05-1.5 µm. The total filler load is 81 %, the total filler volume is 65 %. Radiopacity of Micro Esthetic is 200 % aluminium (radiopacity of enamel: ~200 % aluminium, of dentine:~100 % aluminium).
The composite meets the requirements of: DIN EN ISO 4049, type 1, class 2, group 1.

Indications/Intended use
Restorations of all cavity classes

Performance features
The performance features of the product meet the requirements of the intended use.

Contraindications
Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group
Persons who are treated during a dental procedure.

Intended users
This medical device should only be used by a professionally trained dental practitioner.

Incompatibility with Other Materials
Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Application
Preparation of Tooth and Color Matching
Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.
Ascertain the tooth shade while teeth are still moist and select the appropriate Micro Esthetic shade.
Micro Esthetic is shaded according to VITA®-Shades.

Cavity Preparation
After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. For pulp protection areas close to the pulp should be covered with a thin layer of BisiCAL or a calcium hydroxide material. Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with Best Etch and subsequent application of a light cure bonding agent (e.g. Bond LC) is recommended. Otherwise use a self etching bonding agent (e.g. My E-Bond LC) before application of the Micro Esthetic.

Application and Curing of Micro Esthetic
Apply Micro Esthetic in the selected shade. Micro Esthetic is easy to model. Place it by using a suitable instrument. Matrix strips may be used.
In case of deep and extended cavities best results are obtained with a light cure low able composite e.g. Micro Esthetic flow-viscous or Bisico Bulk Fill, followed by the application of Micro Esthetic.
In case of small cavities Micro Esthetic can be applied directly.
For an optimum result apply Micro Esthetic in layers of max. 2 mm. Light cure each increment separately.

For an optimum result apply Micro Esthetic in layers of max. 2 mm. Light cure each increment separately.
By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm2, cure each increment as follows:

Lighter shades (e.g. A1, A2, C2) 20 seconds

Darker shades (e.g. A3.5, B3) 30 seconds

Hold the light emission window as close as possible to the filling material.

Finishing
Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs. Polish to high gloss with polishing discs or polishing brushes.

Storage
Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Additional Notes/Warnings
- Appli-Tips are for single use only due to hygienic reasons.
- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of the composite restorative.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.
- Keep away from children!

Composition
Dimethacrylates, dental glass, silicon dioxide, photo initiators.

Disposal
Disposal of the product according to local authority regulations.

Reporting obligation
Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note
The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – https://ec.europa.eu/tools/eudamed).

Warranty
Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico GmbH's sole obligation shall be repair or replacement of the Bisico GmbH product.

Limitation of Liability
Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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Bisico GmbH, Johannesswerkstraße 3, 33611 Bielefeld, Germany