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Regidur i®

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Recommendations for use **EN**

Regidur i®

Addition curing extremely hard and fast setting vinyl impression silicone for bite registration

With **Regidur i®** , you have an extremely hard and quick-hardening silicone at your disposal which meets the daily requirements in all respects:

Regidur i® can also be used for sealings of total prosthetic registrations, e.g. when making intraoral pressure-dependent registrations or registrations by means of a face bow.

Regidur i® meets the requirements of ISO 4823, type B (bite registration materials).

Indications/Intended use
- Bite registration material

Performance features
The performance features of the product meet the requirements of the intended use.

Contraindications
None known

Patient target group
Persons who are treated during a dental procedure.

Intended users
This medical device should only be used by a professionally trained dental practitioner.

INSTRUCTIONS FOR USE

Preparing of the cartridge
Remove and dispose the closure-cap and replace it by a supplied mixing cannula. Mount the cartridge into the application gun and extrude needed material.

Bite Registration
Regidur i® should preferably be applied to the patient's lower dental arch. The patient closes their mouth whilst ensuring that the correct occlusion is achieved.

By means of a BRT (Bite-Registration-Tip), Regidur i® can be applied to the closed row of teeth from the outer position to the left or to the right. Particularly in case of abnormal occlusion conditions like for example a cross bite, this feature can be quite advantageous.

After a period of about 45 seconds, the material has reached its final hardness and can be removed from the mouth. Excessive material can be removed with a knife and trimmed with a milling cutter or a stone.

Disinfection
Bite registrations made of Regidur i® can be treated with all disinfectants generally used in the dental area (e.g. 2 % glutaraldehyde solution).

Storage
The material can be stored at room temperature (20–23 °C). Higher storage temperatures reduce the reaction time, lower temperatures prolong it. Maximum storage temperature: < 25 °C

Important notes:
Only very rarely, allergies or other intolerances may occur. In this case, Regidur i® should not be used!

Before first use, insert the new cartridge into the dispenser, remove the sealing cap (do not put the sealing cap onto the cartridge again!) and squeeze out a small quantity of the material until the material comes out of both openings. Only now insert the mixing tip and apply Regidur i®.

The mixing tip used remains on the cartridge until the next application.

Only just before the next application, the mixing tip used is removed and replaced by a new one!

Polysulfide impression materials, ferruginous(III) retraction fibres and solutions as well as constituents of latex gloves may impair the reaction of Regidur i®. To avoid this, the areas concerned should be rinsed thoroughly. (This list makes no claim to be complete.)

Do not use Regidur i® for general dental impression purposes. Its final hardness may cause problems when removing the impression from the mouth.

Do not use Regidur i® after date of expiry!

Please follow the instructions for use of the dispenser!

Disposal
Disposal of the product according to local authority regulations.

Reporting obligation
Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Warranty
Bisico, Bielefelder Dentalsilicone GmbH & Co. KG warrants this product will be free from defects in material and manufacture. Bisico, Bielefelder Dentalsilicone GmbH & Co. KGmakes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico, Bielefelder Dentalsilicone GmbH & Co. KG 's sole obligation shall be repair or replacement of the Bisico, Bielefelder Dentalsilicone GmbH & Co. KG product.

Limitation of Liability
Except where prohibited by law, Bisico, Bielefelder Dentalsilicone GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Technical data
Linear shrinkage rate: max. 0.1 %
Final hardness: approx. 90 Shore-A
Mixing time: not applicable!
Processing time (from beginning of mixing): 15 seconds
Hardening time in the mouth (from beginning of mixing): approx. 45 seconds

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Bisico • Bielefelder Dentalsilicone GmbH & Co. KG • Johanneswerkstraße 3 • D-33611 Bielefeld